

Risk Assessment for School Science Activities

Name and nature of activity	Operating a pressure cooker and autoclave	
Location and date of activity		
Name of teacher/technician	Science ASSIST example risk assessment	
Activity type	<input checked="" type="checkbox"/> Technician procedure <input type="checkbox"/> Teacher demonstration <input type="checkbox"/> Student activity – Student year group	
Physics and general equipment	Type of hazard	Controls and other measures
Electric pressure cooker Autoclave Hotplate (for non-powered steriliser or pressure cooker) Use heat resistant equipment including borrosilicate glassware.	<input type="checkbox"/> Radiation ionising laser <input checked="" type="checkbox"/> Electrical <input checked="" type="checkbox"/> Thermal <input type="checkbox"/> Projectiles <input type="checkbox"/> Sharps <input checked="" type="checkbox"/> Other – Pressure vessel	<input checked="" type="checkbox"/> Relevant signage <input type="checkbox"/> Perspex safety shield <input type="checkbox"/> Sharps container <input checked="" type="checkbox"/> Glassware free from cracks or chips <input checked="" type="checkbox"/> Safety glasses <input checked="" type="checkbox"/> Thermally insulated gloves <input checked="" type="checkbox"/> Other – Ensure that the pressure vessel is not under pressure when opening.
Chemicals used and produced	Type of hazard	Controls and other measures
NA	<input type="checkbox"/> Explosive <input type="checkbox"/> Flammable <input type="checkbox"/> Oxidising <input type="checkbox"/> Gases under pressure <input type="checkbox"/> Corrosive <input type="checkbox"/> Acute toxicity <input type="checkbox"/> Health hazards <input type="checkbox"/> Chronic health hazards <input type="checkbox"/> Environmental <input type="checkbox"/> Other –	<input type="checkbox"/> Limit quantity/concentration <input type="checkbox"/> Perspex safety shield <input type="checkbox"/> Ventilation: natural/exhaust <input type="checkbox"/> Fume cupboard <input type="checkbox"/> Safety glasses <input type="checkbox"/> Laboratory coat/apron <input type="checkbox"/> Gloves: latex/nitrile/neoprene/PVC <input type="checkbox"/> Safety shower <input type="checkbox"/> Other –
Biological/geological materials used	Type of hazard	Controls and other measures
Microbial cultures	<input checked="" type="checkbox"/> Biohazard <input type="checkbox"/> Dust/aerosols <input type="checkbox"/> Sharps <input type="checkbox"/> Manual handling <input type="checkbox"/> Other –	<input checked="" type="checkbox"/> Steriliser <input type="checkbox"/> Disinfectant <input type="checkbox"/> Sharps container <input type="checkbox"/> Dust mask <input type="checkbox"/> Safety glasses <input type="checkbox"/> Gloves <input type="checkbox"/> Other –
Waste produced	Waste disposal procedure	
Sterilised petri dishes containing agar. Excess agar.	<input checked="" type="checkbox"/> Pre-treatment of waste – All microbial and agar waste materials should be sterilised and double bagged prior to disposal in regular waste <input type="checkbox"/> Sink with water – <input checked="" type="checkbox"/> Regular waste – <input type="checkbox"/> Licenced hazardous waste company – <input type="checkbox"/> Other –	
Standard Operating Procedures		
<input checked="" type="checkbox"/> I have read the relevant Standard Operating Procedure. <input checked="" type="checkbox"/> I am experienced/trained in using all the equipment listed. <input checked="" type="checkbox"/> All chemicals used and produced are approved for use. <input checked="" type="checkbox"/> I have read the current SDSs for all hazardous chemicals used and produced. <input checked="" type="checkbox"/> I am aware of safety guidelines for using all chemicals, materials and equipment. <input checked="" type="checkbox"/> I will follow local guidelines for waste disposal (water authority, local council, EPA). <input checked="" type="checkbox"/> I am aware of first aid procedures if required.		
Other comments: Pressure cookers and autoclaves are used for sterilising and do not produce waste as such, but all microbial and agar waste materials should be sterilised and double bagged before disposing in the regular waste.		
Conclusion:		
<input type="checkbox"/> Risks not significant now and not likely to increase. <input checked="" type="checkbox"/> Risks significant but effectively controlled at the moment. <input type="checkbox"/> Risks significant and not adequately controlled at the moment. <input type="checkbox"/> Uncertain about risks; more detailed assessment required.		
Assessment carried out by: Science ASSIST	Signature:	Date: March 2016
Assessment approved by:	Signature:	Date:
Next assessment due:		
This Risk Assessment assumes that the activity will be conducted in a science teaching area with the following facilities: electricity, running water, emergency shut-offs for electricity, gas if applicable, and water, regular testing and tagging of portable appliances; emergency contingencies such as evacuation/emergency plans, appropriate fire extinguishers, spill kits, hand washing facilities, eyewash/safety shower and first aid supplies. It is also assumed that all the necessary licencing requirements and approvals are obtained prior to the activity.		